

# PARAVAN®

MOBILITY FOR LIFE

PARAVAN WHEELCHAIR SERIES

## Order form PR 35 S

Specialist retailer

Commission

HMV no. in 18.99.03.2038



Contact details/customer details

Fill out privacy policy!

**Purchaser/specialist retailer**

**End customer details**

Company

Surname, first name

Contacts

Street/house number

Position

Postcode/city

Street/house number

Telephone/mobile number

Postcode/city

E-mail

Telephone/mobile number

**Advisory sales force**

E-mail

Surname

Additional information

Fill out privacy policy!

**Customer/wheelchair**

**Legal representation**

Medical condition

Surname, first name

Current WC model

Street/house number

**Motor vehicle fitting required**

Yes

No

Postcode/city

Type of fitting

Passive

Active

Telephone/mobile number

**Paying agency**

Current vehicle

Designation

Year of manufacture

Street/house number

Modification of the vehicle so far

Rear entry

Side entry

Postcode/city

Telephone number/e-mail

## Measurements

### Customer measurements

Date of birth

Height

Weight

M1/back height/WC seat

M1/back height/person

M2/seat width/WC seat

M2/seat width/person

M3/seat depth/WC seat

M3/seat depth/person

M4/shoulder width

M5/chest width

M6/seated height in lowest position

M7/distance from seat cushion to upper edge of arm cushion

M8/front edge of arm cushion – measured from backrest

M9/distance from front edge of arm cushion to centre of joystick

M10/lower leg length



M11/foot length

M12/eye measurement in lowest seat position (eye level height)

M13/knee height when seated in wheelchair in lowest position

Measurements measured in customer

Measurements measured in Paravan wheelchair



Notes/drawings/comments



Large empty grey rectangular area for notes, drawings, or comments.

Disclaimer

We, the undersigned, hereby confirm that we have provided the correct measurements in this document and have checked them appropriately. The dimensions are binding for the order.

City, date

Wheelchair user

Signature

Legal guardian

Signature

Enclosures:

Customer photos  
Home visit report  
Data protection declaration  
declaration of assignment




Purchaser/medical supply store/  
company stamp

Medical product advisor

Signature










# Order form | PARAVAN PR 35 S

## Wheelchair – basic equipment

!	Item no.	Picture	Description	Price (excluding VAT)	Selection
	10838	<p><b>PARAVAN PR 35 S</b> Electric wheelchair</p>  <p>Item number for aid: <b>18.99.03.2038</b></p>	<ul style="list-style-type: none"> <li>• Electric standing wheelchair with front wheel drive</li> <li>• Electric standing function</li> <li>• Electric reclining function</li> <li>• Electric lower-leg angle adjustment</li> <li>• Electric back adjustment</li> <li>• Electric 45° tilt backward</li> <li>• Headrest</li> <li>• Seat cushion, contur shape</li> <li>• Seating width 27 - 42 cm</li> <li>• Back cushion, contur shape</li> <li>• Memory function</li> <li>• Very maneuverable, as only 57 cm wide</li> <li>• Ground clearance of around 8 cm</li> <li>• Front tilt supports</li> <li>• Spring-mounted foot rests that can be folded up individually</li> <li>• Calf pads – height and width positions adjustable</li> <li>• Puncture-proof tires with grooved tread of 3 x 8 inches at the front and 200 x 50 at the back with independent spring-mounted wheel suspension</li> <li>• 2 batteries (maintenance-free) - 55 Ah 12V</li> <li>• Pivoting control panel holder</li> <li>• Speed of up to 6 km/h</li> <li>• Transport lock 2 x front and 1 x rear</li> <li>• Maximum user weight: 80 kg</li> </ul>	€ 25,980.00	
	10865	<p><b>Electric standing function</b></p> 	<ul style="list-style-type: none"> <li>• Ergonomic standing</li> <li>• Individually programmable motion sequences</li> </ul>	✓	
	10866	<p><b>Electric reclining function</b></p> 	<ul style="list-style-type: none"> <li>• Sliding back and leg unit with full length compensation</li> </ul>	✓	
	10867	<p><b>Electric back adjustment</b></p>	<ul style="list-style-type: none"> <li>• Biometrically running along with the seat</li> </ul>	✓	

# Order form | PARAVAN PR 35 S

## Wheelchair – basic equipment

!	Item no.	Picture	Description	Price (excluding VAT)	Selection
	10868	<b>Seat tilt</b> 	<ul style="list-style-type: none"> <li>• Seat tilt function</li> <li>• 45° tilt backward</li> </ul>	✓	
	10869	<b>Electric lower-leg angle adjustment</b>	<ul style="list-style-type: none"> <li>• Electric lower-leg angle adjustment</li> <li>• 0-90 degree adjustability</li> <li>• With length compensation</li> </ul>	✓	
	-	<b>Lower-leg length</b>	<ul style="list-style-type: none"> <li>• Lower-leg length                      cm.</li> </ul>	✓	
	61976 61977 16009	<b>Drive options</b>	<ul style="list-style-type: none"> <li>• 6 Km/h</li> <li>• 10 Km/h - incl. TÜV certification (German Road Traffic Licensing Regulations - StVZO), incl. rear sign (<b>mandatory: lighting and seat belt, Gyro control</b>)</li> <li>• 10km/h Export (<b>Mandatory: Gyro control</b>)</li> </ul>	✓ € 880.00 € 340.00	
		<b>Choice of colors</b> Standard			
	10850		• Lemon yellow	✓	
	10851		• Strawberry red	✓	
	10852		• Fairy purple	✓	
	10853		• Unicorn pink	✓	
	10338		• Smurf blue	✓	
	10854		• Pirate blue	✓	
	11985		• Dinosaur green	✓	
	61898	<b>Choice of colors</b> Special color	• Please state RAL/color code	€ 180.00	
		<b>Rims painting</b>			
	13295		• Paintwork in the same colour as cladding	✓	




# Order form | PARAVAN PR 35 S

## Wheelchair – controls / control panels

!	Item no.	Picture	Description	Price (excluding VAT)	Selection
	73055 73062	<b>Control / control panel CJSM</b> 	<ul style="list-style-type: none"> <li>• Right position</li> <li>• Left position</li> </ul>	 ✓ ✓	
	73090 73093	<b>Control / control panel CJSM 2</b> 	<ul style="list-style-type: none"> <li>• Right position</li> <li>• Left position</li> </ul>	 € 780.00 € 780.00	
	73091 73094	<b>Control / control panel CJSM 2 inclusive Bluetooth and infrared</b> 	<ul style="list-style-type: none"> <li>• Position right Bluetooth &amp; Infrarot</li> <li>• Position left Bluetooth &amp; Infrarot</li> </ul>	 € 945.00 € 945.00	
	71860	<b>Control / control panel Gyro control</b>	<ul style="list-style-type: none"> <li>• Electronic track stabilization</li> </ul>	€ 1,090.00	
	--	<b>Control / control panel Special controls</b>	<ul style="list-style-type: none"> <li>• Environmental controls</li> <li>• Special controls</li> <li>• Voice controls</li> </ul>	Upon request	
	--		<ul style="list-style-type: none"> <li>• From front edge of armrest to centre of joystick   STANDARD 100mm</li> </ul>	✓	
	39628	<b>Position Length Joystick Control unit</b>	<ul style="list-style-type: none"> <li>• From front edge of armrest to centre of joystick   80mm</li> </ul>	€ 120.00	
	4823		<ul style="list-style-type: none"> <li>• From front edge of armrest to centre of joystick   120mm</li> </ul>	€ 120.00	

# Order form | PARAVAN PR 35 S









## Wheelchair – controls / control panels

!	Item no.	Picture	Description	Price (excluding VAT)	Selection
	--		<ul style="list-style-type: none"> <li>• Housing at the same height as arm pads   Sleeve STANDARD 60 mm</li> </ul>	✓	
	13783	<b>Height position Joystick control unit</b>	<ul style="list-style-type: none"> <li>• Housing 2 cm above arm pad   Sleeve 40 mm</li> </ul>	€ 170.00	
	13784		<ul style="list-style-type: none"> <li>• Housing 2 cm below arm pad   Sleeve 80 mm</li> </ul>	€ 170.00	
	13785		<ul style="list-style-type: none"> <li>• Housing 4 cm below arm pad   Sleeve 100 mm</li> </ul>	€ 170.00	
	12245		<b>Joystick arrangement for switching sides</b>	<ul style="list-style-type: none"> <li>• Cable left and right</li> <li>• Incl. sleeve for holding the swivel arm</li> </ul>	€ 280.00
	11984	<b>Attendant control R-Net</b> 	<ul style="list-style-type: none"> <li>• Scoot Control with attachment to the backrest</li> </ul>	€ 3,200.00	
	73265	<b>Attendant control R-Net</b> 	<ul style="list-style-type: none"> <li>• Joystick control</li> </ul>	€ 1,260.00	
	35742	<b>Control / control panel R-Net with separate remote control</b>	<ul style="list-style-type: none"> <li>• CJSM</li> </ul>	€ 975.00	
	35743		<ul style="list-style-type: none"> <li>• CJSM 2</li> <li>• STANDARD installation between armrest holder and joystick</li> </ul>	€ 975.00	









# Order form | PARAVAN PR 35 S

## Wheelchair – controls / control panels

!	Item no.	Picture	Description	Price (excluding VAT)	Selection
	61986	<p><b>Tetra fork</b> vertical</p> 	<p><b>Width:</b></p> <p><b>Height:</b></p>	€ 156.00	
	5678	<p><b>Tetragabel</b> horizontal</p> 	<ul style="list-style-type: none"> <li>• with flexible shaft 6 mm (1/4")</li> <li>• 76 mm wide, 83 mm high</li> </ul>	€ 189.00	
	5680	<p><b>Ball attachment</b> foam, red</p> 	<ul style="list-style-type: none"> <li>• large</li> <li>• 68 mm</li> </ul>	€ 76.00	
	27699	<p><b>Ball attachment</b> foam, yellow</p> 	<ul style="list-style-type: none"> <li>• large</li> <li>• 68 mm</li> </ul>	€ 76.00	
	61988	<p><b>Ball attachment</b> foam, blue</p> 	<ul style="list-style-type: none"> <li>• large</li> <li>• 68 mm</li> </ul>	€ 76.00	
	61989	<p><b>Ball attachment</b> plastic</p> 	<ul style="list-style-type: none"> <li>• small</li> <li>• 32 mm</li> </ul>	€ 76.00	
	30181	<p><b>Ball attachment</b> foam, blue striped</p> 	<ul style="list-style-type: none"> <li>• small</li> <li>• 42 mm</li> </ul>	€ 76.00	
	33324	<p><b>Ball attachment</b> Golf ball</p> 		€ 88.00	

## Wheelchair – seat systems / seat part

!	Item no.	Picture	Description	Price (excluding VAT)	Selection
	10871 10872 10873	<b>Back cushion</b> Contoured with lateral support. Choice of sizes (please specify)	<ul style="list-style-type: none"> <li>• Backrest width 27,0 cm - Backrest height 34,0 cm</li> <li>• Backrest width 34,0 cm - Backrest height 38,0 cm</li> <li>• Backrest width 38,0 cm - Backrest height 42,0 cm</li> </ul>	✓	
	46955	<b>Rückenkissen Netti kyphotic</b> 	<ul style="list-style-type: none"> <li>• width 35 cm, height 50 cm</li> </ul>	€ 1,080.00	
	10874 10875 10876	<b>Seat cushion</b> Contour shape Choice of sizes (please specify)	<ul style="list-style-type: none"> <li>• Seat width 27,0 cm - Seat depth 33,0 cm</li> <li>• Seat width 34,0 cm - Seat depth 34,0 cm</li> <li>• Seat width 38,0 cm - Seat depth 42,0 cm</li> </ul>	✓	
	75627	<b>Adjustable side bolsters</b>	<ul style="list-style-type: none"> <li>• on seat left + right</li> <li>• Height and side adjustable</li> </ul>	€ 998.00	
	18830	<b>ADK cushions available</b> (please specify)	<ul style="list-style-type: none"> <li>• Brand / type</li> <li>• Dimensions (WxHxD)</li> </ul>	✓	
	27036	<b>Seat cushion</b> Inko cover 	<ul style="list-style-type: none"> <li>• Inko cover made of Dartex for seat cushion, black</li> </ul>	€ 330.00	
	2113	<b>Calf pad cover</b> Genuine leather	<ul style="list-style-type: none"> <li>• Calf pad entirely in black leather</li> </ul>	€ 252.00	
	13301 13016 13017	<b>Chest strap</b> Monoflex upper body fix for Biolution seat system 	<ul style="list-style-type: none"> <li>• Neoprene chest strap size XS</li> <li>• Neoprene chest strap size S</li> <li>• Neoprene chest strap size M</li> </ul>	✓	
	10884	<b>Kopfstütze</b> Standard 	<ul style="list-style-type: none"> <li>• Height and angle adjustable</li> </ul>	✓	
	71893	<b>Headrest</b> Lateral 	<ul style="list-style-type: none"> <li>• Headrest with lateral side guide</li> <li>• Height and angle adjustable</li> </ul>	€ 399.00	
	62004	<b>Pushing aid</b> with grips 	<ul style="list-style-type: none"> <li>• Pushing aid with grips, mounted on the backrest</li> <li>• Not in combination with holder for lightweight ramp</li> </ul>	€ 280.00	


# Order form | PARAVAN PR 35 S

## Wheelchair – seat systems / seat part / armrests

!	Item no.	Picture	Description	Price (excluding VAT)	Selection
	62005	<b>Pushing aid with grips (customised)</b> 	<ul style="list-style-type: none"> <li>• Pushing aid with grips in the customised design, mounted on the backrest</li> <li>• In combination with warning triangle (10 km/h drive)</li> </ul>	€ 366.00	
	43871	<b>Hooks for adapter plate</b> 	<ul style="list-style-type: none"> <li>• Bag holder</li> </ul>	€ 80.00	
	11789 11790	<b>Armrest holders Variable</b> 	<ul style="list-style-type: none"> <li>• Armrest holder without armrest</li> <li>• Mounted on the backrest as a running unit</li> <li>• 6-axis with 12-fold adjustment, can be swiveled up</li> </ul>	✓ (L) ✓ (R)	
	11787 11788	<b>Armrest Light, short arm shell</b> 	<ul style="list-style-type: none"> <li>• Armrest with light arm shell – 24 cm (short)</li> </ul>	✓ (L) ✓ (R)	
	13302 13303	<b>Armrest Fabric, short</b> 	<ul style="list-style-type: none"> <li>• Padded armrest (fabric) – short 24 cm</li> </ul>	€ 104,00 (L) € 104,00 (R)	
	43303 43304	<b>Armrest Fabric, long</b> 	<ul style="list-style-type: none"> <li>• Padded armrest (fabric) – long 32 cm</li> </ul>	€ 119.00 (L) € 119.00 (R)	
	61685 61684	<b>Armrest Genuine leather, short</b> 	<ul style="list-style-type: none"> <li>• Padded armrest (leather) – short 24 cm</li> </ul>	€ 149.00 (L) € 149.00 (R)	
	61687 61686	<b>Armrest Genuine leather, long</b> 	<ul style="list-style-type: none"> <li>• Padded armrest (leather) – long 32 cm</li> </ul>	€ 169.00 (L) € 169.00 (R)	
	60224 60292	<b>Armrest Light, long arm shell</b> 	<ul style="list-style-type: none"> <li>• Armrest with light arm shell – 32 cm (long)</li> </ul>	€ 76.00 (L) € 76.00 (R)	








# Order form | PARAVAN PR 35 S

## Wheelchair – armrests / arm supports / therapeutic table / footrests

!	Item no.	Picture	Description	Price (excluding VAT)	Selection
	61677 61674	<b>Arm support</b> Synthetic leather, short	<ul style="list-style-type: none"> <li>Cushioned arm support (synthetic leather) – short 24 cm</li> </ul>	€ 104.00 (L) € 104.00 (R)	
	61683 61680	<b>Arm support</b> Synthetic leather, long	<ul style="list-style-type: none"> <li>Cushioned arm support (synthetic leather) – long 32 cm</li> </ul>	€ 119.00 (L) € 119.00 (R)	
	61608 61609	<b>Armrests</b> Reinforcement 	<ul style="list-style-type: none"> <li>Armrest reinforcement – only for long armrests</li> <li>Without edge protection</li> </ul>	€ 56.00 (L) € 56.00 (R)	
	61610 61607	<b>Armrests</b> Extra reinforcement 	<ul style="list-style-type: none"> <li>Armrest reinforcement – only for long armrests</li> <li>With edge protection</li> </ul>	€ 112.00 (L) € 112.00 (R)	
	61613 61611	<b>Armrest holders</b> 8 cm extension 	<ul style="list-style-type: none"> <li>Armrest extension by 8 cm on left side</li> <li>Armrest extension by 8 cm on right side</li> </ul>	€ 88.00 (L) € 88.00 (R)	
	61614 61612	<b>Armrest holders</b> 12 cm extension 	<ul style="list-style-type: none"> <li>Armrest extension by 12 cm on left side</li> <li>Armrest extension by 12 cm on right side</li> </ul>	€ 108.00 (L) € 108.00 (R)	
	69496	<b>Therapeutic table</b> 	<ul style="list-style-type: none"> <li>Therapeutic table, installation on armrest</li> <li>Foldable and removable</li> <li>Always mount on the opposite side of the joystick positioning</li> </ul>	€ 890.00	
	69494	<b>Therapeutic table</b> Barrier	<ul style="list-style-type: none"> <li>Table with barrier (to prevent objects from slipping down)</li> <li>PLEASE NOTE: Please always fill in the additional "Dimension sheet for therapeutic table" form</li> </ul>	€ 970.00	
	10858	<b>Footrests</b> manually foldable, individual 	<ul style="list-style-type: none"> <li>Footrests – manually foldable, individual, spring-loaded</li> <li>Footrest with a length of 20 cm</li> </ul>	✓ narrow	





# Order form | PARAVAN PR 35S

## Wheelchair – wheelchair safety / lighting / tires

!	Item no.	Picture	Description	Price (excluding VAT)	Selection
	44509	<b>Locking plate Dahl docking station</b> 	<ul style="list-style-type: none"> <li>• Plate incl. installation</li> </ul>	€ 820,00	
	62020	<b>Headlights</b> 	<ul style="list-style-type: none"> <li>• Frontal headlight (black)</li> <li>• Rear light emitting diodes, including turn signals</li> </ul>	✓	
	73195	<b>Tyres</b> Grooved tread PR 	<ul style="list-style-type: none"> <li>• Grooved tread, puncture-resistant, black</li> <li>• front and rear axle</li> </ul>	✓	
	73170	<b>Tyres</b> Lug tread PR 	<ul style="list-style-type: none"> <li>• Lug tread, puncture-resistant, black</li> <li>• Front axle (both tyres)</li> </ul>	€ 45.00	
	13049	<b>Tyres</b> Lug tread PR 	<ul style="list-style-type: none"> <li>• Lug tread, puncture-resistant, black</li> <li>• Rear axle (both tyres)</li> </ul>	€ 50.00	
	73211	<b>Tyres</b> Grooved tread pneumatic 	<ul style="list-style-type: none"> <li>• Grooved tread, pneumatic tyres, black</li> <li>• Front axle (both tyres)</li> </ul>	€ 69.00	
	73210	<b>Tyres</b> Lug tread pneumatic 	<ul style="list-style-type: none"> <li>• Lug tread, pneumatic tyres, black</li> <li>• Front axle (both tyres)</li> </ul>	€ 69.00	

# Order form | PARAVAN PR 35 S

## Wheelchair – general accessories / batteries and charging

!	Item no.	Picture	Description	Price (excluding VAT)	Selection
	62023	<b>Power supply</b> 24 volt	<ul style="list-style-type: none"> <li>• Cable around 50 cm for further use</li> </ul>	€ 238.00	
	41069	<b>Double USB charging socket</b> 	<ul style="list-style-type: none"> <li>• External charging socket for charging</li> <li>• Please specify position:</li> </ul>	€ 569.00	
	71955	<b>Disconnecters</b> Battery	<ul style="list-style-type: none"> <li>• Battery disconnecter for transport in planes, in the middle of the chassis at the back</li> </ul>	€ 347.00	
	62006 62007	<b>Small bag</b> 	<ul style="list-style-type: none"> <li>• Including installation on left armrest</li> <li>• Including installation on right armrest</li> </ul>	€ 156.00 € 156.00	
	12317	<b>Backpack</b> 	<ul style="list-style-type: none"> <li>• Paravan backpack</li> </ul>	€ 198.00	
	12967 12968 12969	<b>Charger</b> 	<ul style="list-style-type: none"> <li>• Charger, without fan DE</li> <li>• Charger, without fan US</li> <li>• Charger, without fan GB</li> </ul>	✓	



Customised design form (addendum to quotation: )

Customised design  
text

Material costs

Sale price/€

Working time/hours

Discount:

Yes  
No

Date/signature of service manager


Production description

Technically checked and released

Date/signature



**PARAVAN therapy table dimensions**



<b>A</b>	Total width		mm
<b>B</b>	Belly section width		mm
<b>C</b>	Belly section depth		mm
<b>D</b>	Depth without belly section		mm
<b>E</b>	Distance to cut-out for centre control		mm (Min. 30 mm)
<b>F</b>	Cut-out for control unit on armrest	Right	Left
	Mounting side Bracket	Right	Left
	Draw in positioning of picobuttons	<b>Please tick in the picture above</b>	

Date / Signature



MOBILITY FOR LIFE

**PARAVAN GmbH**

Paravanstraße 5-10

72539 Pfronstetten-Aichelau

Germany

Tel.: +49 (0) 73 88 / 99 95-612

Fax: +49 (0) 73 88 / 99 95-999

E-mail: [info@paravan.com](mailto:info@paravan.com)