



Crash tested  
according to the  
latest guidelines



# DOCKING- STATION

The perfect connection between  
car and wheelchair

Version Dockingstation V5.0 09/2020 EN

- Top production and product quality, "Made in Germany"
- Mechanical or electrical locking system
- Emergency release device
- Can also be approached at an inclined angle due to predetermined guide
- 100% safety and stability (stable in all three axes)
- No instability is transmitted and it is without momentum, therefore driving with SPACE DRIVE® is safe and smooth
- High strength (electro-galvanised) steel is used to prevent any torsion
- For self-drivers and passengers
- Developed entirely at PARAVAN
- Manufactured entirely at PARAVAN
- Can be installed in many vehicles
- Ground clearance under wheelchair approx. 75 mm
- Crash-tested in compliance with the most recent guidelines at 20G
- Crash-tested wheelchair locking set
- 24-hour spare parts delivery