

# The Paravan PR 25 wheelchair



Read in the QR code and go directly to the product video: 

## *As flexible as no other powerchair*

With our newest member of the PARAVAN power wheelchair series, you can comfortably lower the seat all the way down to the ground. You can easily reach objects on the floor and children can play with others and be at the same height for story or play time. Another

unique feature is the ability to move the seat completely forward horizontally, which makes transfers a lot easier. In addition, the PR 25 is equipped as standard with an easy-to-use joystick control and an individually adjustable seat system that grows with the child.



Photos may contain optional extras.

## Your benefits:

### Series model PR 25:

- Power wheelchair with seat-to-floor function
- Entry height of only 15 cm when seat is on ground
- Driving in lowered position possible
- Driving when seat is forward possible
- Stepless seat adjustment in length and lowering
- Highly manoeuvrable thanks to width of only 64cm
- Adaptive PARAVAN seat shell concept
- R-Net CJSM joystick control
- Puncture proof tyres with grooved tread of 3 x 8 inches at the front and 200 x 50 inches with independent sprung suspension
- 10 A outdoor charger without fan
- Swing away control panel holder
- Seat and backrest width extension
- Standard speed 4 mph, 6.5 mph optional
- Transport tie down points 2 x front and 1 x rear
- Max. user weight: 60 kg
- Two-way spring-mounted chassis at the rear with front tilt support
- 12-way adjustable and rotating armrest mounts
- 85 AH batteries
- Black tyres (non-marking)

### Key facts and options for PR 25

- Made in Germany: Highest manufacturing and product quality developed entirely by Paravan
- 24-hour spare part service
- Bespoke assembly possible to meet individual needs
- Modular design which can be adapted to meet the needs of the individual clinical requirements
- Large choice of individual seat pads (optional)
- Gyro tracking control (optional)
- Optional controls available
- Special controls and environmental controls available
- Individual, modern paintwork options
- Vehicle docking to secure the wheelchair in a vehicle (optional)
- Electric tilt angle up to 45° backwards and 20° forwards (optional)
- Seat lift (optional)
- 3- or 4-point seat belt possible



#### Highlight: Seat-to-floor function

Thanks to the electric length and lowering function, the seat can move forward and downward all the way to the floor.

#### Highlight: Getting in and out from the floor

With the PR 25 we offer a unique power chair with an entry height of just 15 cm when the seat is on the floor.