

PARAVAN PR 35 S

Only the best for the little ones

Our PR 35 S is a very special power wheelchair. It is small, compact, colourful and versatile. In other words, ideally developed and tailored to the needs of our smallest customers. Lying down, standing up, tilting - all functions are on board. And the 30 percent smaller chassis also makes the wheelchair extremely agile. Just the way children love it! Whether at day care, school or at home - our PARAVAN PR 35 S is made to be actively involved. Right in the thick of the action!



YOUR BENEFITS



Highlight: Standing, sitting, lying & tilting

Modern, flexible power wheelchair with comprehensive features as standard for children or smaller adults.

Highlight: Growing seat concept

The seat concept of the PARAVAN PR 35 S grows with the user, making the wheelchair the perfect companion from childhood to teenage years.

PR 35 S basic model:

- Power wheelchair with front wheel drive
- Integrated standing, lying and tilting function.
 Standing position directly from sitting or lying position.
- Very manoeuvrable, as only 570 mm wide
- Adaptive PARAVAN seat shell concept
- Joystick control R-Net with memory function
- Electric lower leg angle adjustment with length compensation
- Electric backrest adjustment, biometrically running along with the seat
- Electric tilt 45° to the rear, 20° to the front
- Footrest spring-loaded and foldable
- Calf position adjustable in height and width
- Puncture-proof tyres with grooved tread 3.00 x 8 inch at the front and 200 x 50 mm at the rear with independent wheel suspension sprung
- Charger Outdoor-10 Amp without fan
- Swing-away control panel holder
- Speed max. 6 km/h (10 km/h optional)
- Transport lock 2x front and 1x rear
- Max. user weight: 80 kg
- 12-way adjustable armrest
- · Batteries 55 AH

Keyfacts and options PR 35S:

- Made in Germany: Best manufacturing and product quality. Completely developed inhouse by PARAVAN
- Custom modifications possible
- Seat system that grows with the user in terms of seat width, seat depth and backrest
- Extremely quiet servomotors ideal for schools or at day care
- Modular construction system, individually adaptable to the needs of the clinical picture
- Large selection of individual seat cushions (optional)
- Electric track stabilisation (optional)
- Special controls and environment controls possible (optional)
- Individual modern paintwork
- Floor fixation for securing in the vehicle (optional)