

PARAVAN PR40

Versatility reinvented

Sometimes the time is ripe to completely rethink things and explore new paths. We have done just this with the development of our new variable and flexible PR 40 power wheelchair – a perfect symbiosis of our previous PR series models. The PR 40 combines the technical perfection of the PR 50 with the multifunctionality of the PR Biolution at unbeatable value for money. As standard, the PR 40 comes with an electric tilt, standing, and lying position and a redesigned adaptive seat shell system for the perfect level of versatility.



YOUR BENEFITS



Highlight: value for money

Modern variable electric standing wheelchair with extensive standard fittings at a fair price

Highlight: delivered within 14 days at the latest

With the PR 40, we are offering a first-class electric wheelchair with guaranteed delivery within 14 days of placing the order (without special assembly)

Series model PR 40:

- Electric standing wheelchair with front wheel drive
- Integrated standing, lying, and tilting function
- Very maneuverable, only 640 mm wide
- Adaptive PARAVAN seat shell concept
- R-Net CJSM joystick control
- Electric lower-leg angle adjustment with length compensation
- · Electric back adjustment, biometrical
- Electric 45° tilt backwards and 20° tilt forwards
- Separate spring-mounted foot rests that can be folded up individually
- Calf pads height and width positions adjustable
- Puncture-proof tires with grooved tread of 3 x 8 inches at the front and 200 x 50 inches with independent suspension sprung
- 10 A outdoor charger without fan
- Swing away control panel holder
- Standard speed 4 mph
- Transport tie down points 2 x front and 1 x rear
- Maximum user weight: 140 kg (weight of person)
- Two-way spring-mounted chassis at the rear with front tilt supports
- 12-way adjustable and rotating armrest
- 85 AH batteries

Key facts and options for PR 40:

- Made in Germany: the highest level of manufacturing and product quality. Developed by PARAVAN entirely in-house
- Special assembly to meet individual needs
- Seat and backrest width extension (optional)
- Modular construction system, can be adapted individually to the disability profile
- Large selection of individual seat cushions (optional)
- Electric track gyro stabilisation (optional)
- Optional further controls available
- Special controls and environmental controls available
- Customised modern paint finishes
- Vehicle docking station to secure the wheelchair in a vehicle (optional)